

MEDISAFE[®] FINETOUCH[®] II


Lancing device for
capillary blood sampling

The
puncture
with a
gentle
touch



 **TERUMO**

Less pain¹



- 5 adjustable penetration depths. Adaptable to individual patient's skin type. Minimum penetration depth is 0.3 mm.
- Hollow needle specially cut for noticeable pain reduction.¹

Protected use



The lancet needle is surrounded by a plastic holder, before and after penetration, preventing needle stick injury.

Easy Handling



The lancet is inserted into the lancing device in one simple move.




The slim lancet holder tip makes the desired puncture site particularly easy to target.




By sliding the ejector the used lancet can be disposed.

Specifications

	Product Name	Code	Pieces
	MEDISAFE FINETOUGH II Lancing Device with Set Lever	MS*FT2R	1 / 10 / 100
	TERUMO FINETOUGH Lancet	Country related [△]	25 / 250 / 2500 [△] 30 / 300 / 3000 [△]

Contains neither components made of natural rubber latex nor PVC. [△] Please contact your local Terumo office.



¹ Intra company data shows that puncture with MEDISAFE FINETOUGH II lancing device is experienced less painful compared to TERUMO FINETOUGH lancing device.

This document is an overview and not intended to be comprehensive. For more details or in case of questions we refer to the instructions for use delivered with the product or call your local Terumo office. MEDISAFE FINETOUGH II lancing device is manufactured by TaiDoc Technology Corporation. FINETOUGH lancet is manufactured by Terumo Corporation. TERUMO®, MEDISAFE® and FINETOUGH® are registered trademarks of Terumo Corporation.



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

EMEA Sales Offices

**Terumo Europe NV
Africa Business Division**
+32 16 38 13 08

**Terumo Europe NV
Benelux Sales Division**
Belgium:
0800 14468
The Netherlands:
0800 0231938

**Terumo Europe NV
Emerging Market Division**
+32 16 38 12 11

Terumo Deutschland GmbH
+49 6196 80 230

**Terumo Deutschland GmbH
Switzerland**
+41 56 419 10 10

Terumo Europe España SL
+34 902 10 12 98

Terumo France S.A.S.
+33 130 96 13 00

Terumo Italia S.r.l.
+39 06 94 80 28 00

Terumo Russia LLC
+7 495 988 4740

Terumo Sweden AB
+46 31 748 58 80

**Terumo Sweden AB
Denmark**
+45 7020 93 80

Terumo Middle East FZE
+971 4 292 0200

Terumo UK Ltd
+44 1276 480 440

**Terumo BCT Tıbbi Cihazlar
Dağıtım ve Hizmetleri A.Ş.**
+90 216 645 92 00



© Registered Trademark
Published by Terumo Europe NV
DI 143GB-0318FK-I(03.18)E